Description

ExcITE® F is a light-curing, nano-filled, fluoride-releasing, single-component adhesive for dentin and enamel bonding in conjunction with the total-etch technique.

Composition

ExcITE F contains phosphonic acid acrylate, HEMA, dimethacrylate, highly dispersed silicone dioxide, initiators, stabilizers and potassium fluoride in an alcohol solution.

Indication

- Adhesive for direct light-curing and dual-curing composite and compomer restorations
- Adhesive for indirect all-ceramic and composite restorations (inlays, onlays, veneers; only light-curing)

Contraindication

- Do not use ExcITE F if a patient is known to be allergic to any of the material's ingredients or if the stipulated working technique cannot be employed.
- ExcITE F must not be used in combination with purely self-curing composite materials. Dual-curing materials always have to be light-activated.
- Do not use ExcITE F for direct pulp capping.

Side effects

In rare cases, ingredients of ExcITE F may cause sensitization. In such cases, the product should cease to be used.

Interactions

Materials containing eugenol or clove oil may inhibit the polymerization of ExcITE F.

Application**Direct restorations**

1. A dry operating field must be ensured. Ideally, a rubber dam should be placed (e.g. OptraDam® Plus).
2. In very deep cavities, areas close to the pulp should be selectively coated with a calcium hydroxide liner (e.g. ApexCal®) and subsequently covered with a pressure-resistant cement (glass ionomer cement, e.g. Vivaglass® Liner).
3. Apply phosphoric acid gel (e.g. Total Etch) to the prepared enamel and then flow the etchant onto the prepared dentin. The etchant should be left to react on the enamel for 15-30 seconds and on the dentin for 10–15 seconds. Following this, remove all etchant gel with a vigorous water spray for at least 5 seconds. Excess moisture should be removed leaving the dentin surface with a glossy wet appearance (wet bonding). This can be done with either an air gun, a dry brush, a foam pellet or other lint-free absorbents. **Do not overdry the dentin!**
4. **Use of the different delivery forms**

ExcITE F is available in bottles, the VivaPen® and Soft Touch™ Single-Dose vessels.

VivaPen delivery form

- Remove the cap of the VivaPen by turning it beyond the point of resistance and simultaneously pulling it off.
- Snap on the VivaPen brush cannula. The brush cannula can be rotated up to approx. 270 degrees to the desired position.

Attention: Do not rotate the cannula beyond the point of resistance; otherwise it will detach!

- Activate the click mechanism several times to saturate the VivaPen brush cannula with ExcITE F. The yellow colour of ExcITE F becomes visible once the brush cannula is saturated.

Soft Touch Single-Dose vessel with pre-mounted applicator brush

Activate the Soft Touch Single-Dose vessel to prepare it for use as follows:

- Place the Single-Dose vessel on a firm surface.
 - Press the pre-mounted applicator brush into the Single-Dose vessel until a noticeable resistance is overcome.
5. Saturate enamel and dentin with a generous amount of ExcITE F and agitate the adhesive on the prepared surfaces for at least 10 seconds. Make sure that all the cavity walls are completely covered.
 6. Disperse ExcITE F to a thin layer with a weak stream of air, thereby removing any excess. Pooling should be avoided. A uniform, glossy appearance shows that the tooth surface is completely sealed.
 7. Light-cure ExcITE F for 10 seconds at a light intensity of more than 500 mW/cm² (with e.g. bluephase® using the LOP mode).
 8. Apply the restorative material (e.g. Tetric EvoFlow® or Tetric EvoCeram®) according to the instructions of the manufacturer.
 9. VivaPen: Discard the used VivaPen brush cannula after disinfecting the VivaPen with disinfectant wipes and replace it with the original cap. The VivaPen brush cannula is removed by turning it beyond the point of resistance and simultaneously pulling it off. A new VivaPen brush cannula must be attached for every application.

Indirect restorations

1. Prepare the surfaces of the composite or ceramic restoration according to the instructions of the manufacturer.
2. Prepare the tooth surfaces according to the procedure indicated for direct restorations, step 1–7.

Important note: Avoid pooling, as this can compromise the accuracy of fit of the definitive restoration.

3. Following this, seat the restoration using a light-/ dual-curing luting material (e.g. Variolink® II) and light-cure according to the instruction of the manufacturer.

Note

- After each use, the VivaPen must be disinfected using disinfectant wipes for hygienic reasons. Additionally, the VivaPen can be protected against contamination with the VivaPen sleeve.
- Do not use oxidizing disinfectants.
- Do not immerse the VivaPen in disinfectant solutions.
- Do not use other types of brush cannulas.
- The fill-level indicator of the VivaPen identifies the remaining amount of adhesive: If the white indicator is no longer visible on the display, the VivaPen is empty.

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Dental adhesive with
fluoride release

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Warnings

- ExcITE F is an irritant. Contact with the skin, mucous membrane or eyes should be avoided. In case of skin contact, immediately wash with copious amounts of water.
- If the material should come into contact with the eyes, immediately rinse with copious amounts of water for at least 15 minutes and consult a physician or ophthalmologist.
- In rare cases, contact with the skin may cause sensitization to the ingredients.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylates.

Note on Total Etch

Total Etch is a gel for enamel etching and dentin conditioning. It contains phosphoric acid (37 wt% in water), thickening agent and colour pigments. The surfaces of glass ionomer-based bases/liners are dissolved by Total Etch.

Warning

-   When filling the Total Etch small syringe always wear protective goggles and gloves. Observe the Total Etch Instructions for Use and do not fill the Total Etch syringe with etchant gel from jumbo syringes of other manufacturers.
-  Total Etch contains phosphoric acid and is corrosive. Avoid contact with the eyes, mucous membrane and skin (goggles are recommended for both patients and operators). After contact with the eyes, immediately rinse them with copious amounts of water for at least 15 minutes and see a physician/ophthalmologist.
- After skin contact, immediately rinse with copious amounts of water.
- After contact with sensitive materials (surfaces, clothing), immediately clean with copious amounts of water.

Storage

- Expiry date: see label/imprint on VivaPen/packaging
- Close bottles/VivaPen with the original cap immediately after use
- Storage temperature: 2–28 °C (36–82 °F)

Keep material out of children's reach.

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

The Material Safety Data Sheet is available online at www.ivoclarvivadent.com

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